

# Silencing Gene Expression: *siRNA Delivery by Transfection*

## Introduction

Delivery of exogenous small interfering RNA (siRNA) in cultured cells is an effective method to modulate target gene expression via RNA interference (RNAi). However, genetic manipulations in primary neuronal cultures are especially inefficient and often toxic. In fact, neurons are considered one of the most difficult and resistant cell types for introduction of siRNA oligonucleotides.

Several siRNA delivery systems and reagents are available depending on the cell type, cell culture preparation, and desired level of target gene silencing. The protocol presented here demonstrates siRNA delivery by transfection for efficient GAPDH gene silencing in iCell® GABANeurons (formerly known as iCell Neurons) using the Accell siRNA technology in 96-well cell culture plates. This Application Protocol serves as a guide for delivery of other Accell siRNAs for use in different plate formats and endpoint readouts.

## Required Equipment and Consumables

The following equipment and consumables are required in addition to the materials specified in the iCell GABANeurons User's Guide.

| Item   | Vendor                                | Catalog Number(s)   |
|--|---------------------------------------|---------------------|
| <b>Equipment</b>   |                                       |                     |
| Multichannel Pipettors and Sterile Tips  | Multiple Vendors                      |                     |
| Quantitative Real-time PCR (qRT-PCR) Instrument  | Multiple Vendors                      |                     |
| <b>Consumables</b>   |                                       |                     |
| iCell GABANeurons Kit, 01279 <sup>1</sup>  | Cellular Dynamics International (CDI) | R1011, R1084, R1118 |
| iCell GABANeurons Kit, 01434 <sup>1, 2</sup>   | Cellular Dynamics International (CDI) | R1013, R1053        |
| 1.5 ml RNase-free Tubes  | Multiple Vendors                      |                     |
| 5X siRNA Buffer  | Thermo Fisher Scientific              | B-002000-UB-100     |
| 96-well Flat-bottom Microplate, TC-treated, Falcon (96-well Cell Culture Plate) <sup>3</sup> | STEMCELL Technologies                 | 38022               |
| Accell siRNA Delivery Media  | Thermo Fisher Scientific              | B-005000-500        |

| Item                                  | Vendor                   | Catalog Number(s) |
|---------------------------------------|--------------------------|-------------------|
| Accell siRNA Oligonucleotides (siRNA) | Thermo Fisher Scientific |                   |
| Nuclease-free Water                   | Multiple Vendors         |                   |
| RNAse-free Water                      | Multiple Vendors         |                   |

1 Order the kit whose iCell GABANeurons were derived from the desired donor. iCell GABANeurons, 01279 and iCell GABANeurons, 01434 were derived from apparently healthy, normal donors. iCell GABANeurons, 01434 are exclusive to CDI.

**Note:** This Application Protocol was optimized using iCell GABANeurons Kit, 01434 (Cat. No. R1013). CDI anticipates you will achieve similar results using iCell GABANeurons derived from other apparently healthy, normal donors.

2 Formerly known as iCell Neurons (Cat. No. NRC-100-010-001).

3 Similar products are available from multiple vendors.

## Workflow

iCell GABANeurons are thawed and plated into a 96-well cell culture plate previously coated with poly-L-ornithine and laminin solutions. On day 4 post-plating, medium is replaced, and cells are transfected. The siRNA-mediated gene silencing is analyzed subsequently at the optimal post-transfection timepoint for the target gene and endpoint assay.



## Methods

### Culturing iCell GABANeurons

1. Coat a 96-well cell culture plate with a base layer of 0.01% poly-L-ornithine solution and a top coating of a 3.3 µg/ml laminin solution according to the iCell GABANeurons User's Guide.
2. Prepare the Complete Maintenance Medium according to the iCell GABANeurons User's Guide.
3. Thaw the neurons according to their User's Guide.
4. Remove a sample of the cell suspension and count the viable neurons using a hemocytometer.
5. Further dilute the cell suspension in Complete Maintenance Medium to 400,000 cells/ml.
6. Aspirate the laminin solution from the 96-well cell culture plate. Immediately add 100 µl/well of the cell suspension (40,000 cells/well).
7. Incubate in a cell culture incubator at 37°C, 5% CO<sub>2</sub>. Maintain the neurons according to their User's Guide until ready to perform the transfection.

### Transfecting iCell GABANeurons with siRNA

The following procedure details transfection of the neurons cultured on day 4 post-plating. The volumes indicated are per well of a 96-well cell culture plate. Scale volumes appropriately for multiple wells or for other well formats.


1. Warm the Accell siRNA Delivery Media to room temperature.
2. Dilute 5X siRNA Buffer in RNase-free water to achieve a 1X siRNA Buffer.
3. Briefly centrifuge each tube containing siRNA to ensure the pellet is at the bottom of the tube.
4. Reconstitute the siRNA in 1X siRNA Buffer to achieve a 100  $\mu$ M stock solution.
5. Dispense 125  $\mu$ l of Accell siRNA Delivery Media in a 1.5 ml RNase-free tube.
6. Add 1.25  $\mu$ l of 100  $\mu$ M siRNA stock solution to achieve a 1  $\mu$ M siRNA delivery mix. Mix gently.

**Note:** Store the remaining 100  $\mu$ M siRNA stock solution at  $-20^{\circ}\text{C}$ .

7. Remove the 96-well cell culture plate containing the neurons from the incubator.
8. Aspirate the Complete Maintenance Medium from the 96-well cell culture plate and replace with 100  $\mu$ l/well of 1  $\mu$ M siRNA delivery mix.
9. Incubate in a cell culture incubator at  $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$  for 24 hours.
10. Aspirate 1  $\mu$ M siRNA delivery mix and replace with 100  $\mu$ l/well of fresh Complete Maintenance Medium.
11. Incubate in a cell culture incubator at  $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$ . Maintain the neurons according to their User's Guide until the optimal timepoint to analyze the siRNA-mediated gene silencing is attained.

### Measuring siRNA-mediated Gene Silencing in iCell GABANeurons

The optimal post-transfection timepoint at which to measure the RNAi-mediated effect depends on the target gene, siRNA, and endpoint assay. In a representative experiment, GAPDH mRNA silencing in iCell GABANeurons (Figure 1) was measured on day 7 post-plating (day 3 post-transfection) using the TaqMan Gene Expression Cells-to- $\text{C}_T$  Kit (Life Technologies, Cat. No. 4399002M) according to the manufacturer's instructions.



**Figure 1: RNAi-mediated Gene Silencing in iCell GABANeurons**

On day 4 post-plating, iCell GABANeurons, 01434 were transfected with Accell siRNA Delivery Media only (mock), a control (scrambled) siRNA in Accell siRNA Delivery Media, or a titration of GAPDH siRNA in Accell siRNA Delivery Media. On day 7 post-plating (day 3 post-transfection), GAPDH mRNA levels were measured relative to 18s rRNA levels and compared to the mRNA levels obtained following transfection of the control siRNA (mean  $\pm$  SEM,  $n = 3$  independent siRNA delivery mixes).

---

## Summary

iCell GABANeurons provide an in vitro test system that recapitulates native human neurobiology while the Accell siRNA technology offers an efficient method for siRNA transfection experiments. The methods presented here highlight the specificity with which specific gene function can be examined precisely in human neurons via RNAi. Furthermore, the data suggest that the impact of RNAi-mediated knockdown of specific target genes will be effective when monitored through other endpoint assays.

**Customer's Responsibilities**

CDI does not guarantee that you will obtain equivalent results from using iCell or MyCell products as described herein or that such use will not infringe any intellectual property right(s) of any third party(ies). You are solely responsible for obtaining any licenses you may require for your specific research use(s) of the iCell or MyCell products not expressly conveyed under CDI's terms and conditions of sale or other transfer of the iCell or MyCell products to you.

**Conditions of Use**

For life science research use only.

**Trademarks**

iCell and MyCell are registered trademarks, and Cellular Dynamics and the  logo are trademarks of Cellular Dynamics International, Inc.

All other brands, product names, company names, trademarks, and service marks are the properties of their respective owners.

**Copyright Notice**

© 2018 Cellular Dynamics International, Inc. All rights reserved.

**Revision History**

Version 1.1: April 2018  
AP-GABAsiRNA180418