

Measuring Cardiac Activity: *Functional Maturation on the xCELLigence RTCA ePacer*

Introduction

iCell Cardiomyocytes² have been optimized for rapid recovery from cryopreservation. Their derivation from human induced pluripotent stem cells, high purity, functional relevance and ease of use combine to make iCell Cardiomyocytes² an optimal in vitro test system for interrogating cardiac biology in basic research and many areas of drug development.

The xCELLigence RTCA ePacer is a scalable system that allows for functional maturation of iCell Cardiomyocytes² by providing electrical pacing with less stress than currently available models. By controlling the electrical activity, organized sarcomeric structure improves and the cells demonstrate a positive force frequency relationship which increases the accuracy of measurements, well to well consistency and experimental alliance. The system provides a high-throughput assay platform compatible with many options of relevant readouts for safety/toxicology assessment and drug discovery.

This Application Protocol describes how to handle and perform pacing using the iCell Cardiomyocytes² on the xCELLigence RTCA ePacer for functional maturation.

Required Equipment, Consumables and Software

The following equipment, consumables and software are required in addition to the materials specified in the iCell Cardiomyocytes² User's Guide.

Item	Vendor	Catalog Number
Equipment		
Multichannel Pipettor, 20 and 200 µl	Multiple Vendors	
xCELLigence RTCA CardioECR	ACEA Biosciences, A part of Agilent Technologies	
xCELLigence RTCA ePacer	ACEA Biosciences, A part of Agilent Technologies	
Consumables		
iCell Cardiomyocytes ² Kit, 01434	FUJIFILM Cellular Dynamics, Inc.	R1017 CMC-100-012-000.5
Centrifuge Tubes, 50ml Sterile	Multiple Vendors	
Dulbecco's Phosphate Buffered Saline without Ca ²⁺ and Mg ²⁺ (DPBS)	ThermoFisher Scientific	14190
E-Plate CardioECR 48	ACEA Biosciences, A part of Agilent Technologies	00300600940 00300601110
Fibronectin	Roche	11051407001 11080938001
Microcentrifuge Tubes	Multiple Vendors	
Sterile Reagent Reservoirs	Multiple Vendors	
Sterile Water	Multiple Vendors	


Software	
RTCA CardioECR Data Acquisition and Analysis Software	ACEA Biosciences, A part of Agilent Technologies
RTCA ePacer Instrument Software	ACEA Biosciences, A part of Agilent Technologies

Workflow

The cardiomyocytes are thawed and plated into an E-Plate previously coated with fibronectin. Plating Medium is replaced with Maintenance Medium 4 hours after plating and Maintenance Medium is replaced every 48 hours thereafter.

Electrical stimulation (pacing) of the cardiomyocytes begins on day 7 at 1.0 Hz and increases by 0.5 Hz every 5 days. During pacing, Maintenance Medium is replaced every 24 hours.

On day 22, force-frequency relationship is measured. Cells can be treated with compounds and cardiac activity recorded when the force-frequency relationship is positive.



Methods

Note: Prior to beginning, place the RTCA CardioECR and RTCA ePacer stations in a cell culture incubator at 37°C, 5% CO₂.

Preparing the E-Plate - Day 1

1. Reconstitute fibronectin in sterile water at 1 mg/ml according to the manufacturer's instructions.
2. Dilute 1 mg/ml fibronectin solution in DPBS to a final concentration of 10 µg/ml.
3. Add 50 µl/well of the 10 µg/ml fibronectin solution to the center of the wells of the E-Plate to evenly coat the bottom of the well.
4. Incubate at 37°C for at least 1 hour.

Thawing Cardiomyocytes- Day 1

1. Equilibrate iCell Cardiomyocytes Plating Medium to 37°C.
2. Aspirate the fibronectin solution from the E-Plate.

Note: Do not allow the fibronectin-coated surface to dry.

3. Add 50 µl/well of Plating Medium to the center of the wells.
4. Equilibrate the E-Plate in a cell culture incubator at 37°C, 5% CO₂ for at least 15 minutes.
5. Transfer the E-Plate to the RTCA CardioECR station. Record a background measurement according to the RTCA CardioECR Instrument Operator's Guide.
6. Insert the E-Plate into RTCA ePacer station. Record a background measurement according to the RTCA ePacer Instrument Operator's Guide.
7. Thaw the cardiomyocytes according to their User's Guide to a final volume of 4 ml. Dilute the 1 ml cell suspension from the cryovial in 1 ml of Plating Medium, rinse, and add an additional 2 ml of Plating Medium.
8. Use the viable cells/vial from the Certificate of Analysis to calculate the final volume of Plating Medium needed to obtain a final cell plating density of 1 x 10⁶ viable cells/ml.
9. Remove a sample of cells to confirm viability using a hemocytometer (using trypan blue dye exclusion). If the cell viability is below 50%, contact Technical Support for assistance.
10. Dispense 50 µl/well of the cell suspension (50,000 cells/well) to the center of each well of the E-Plate.
11. Leave the E-Plate undisturbed at room temperature for approximately 30 minutes to allow the cells to settle.
12. Culture the cardiomyocytes in a cell culture incubator at 37°C, 5% CO₂ for 4 hours.
13. Equilibrate an aliquot of Maintenance Medium to 37°C.
14. Replace the Plating Medium with Maintenance Medium. Tilt the E-Plate, remove the spent medium, and gently add 100 µl/well of Maintenance Medium to the side of the well to avoid disturbing the cells.
15. Culture the cardiomyocytes in the RTCA ePacer station at 37°C, 5% CO₂.

Note: Inserting the E-Plate into the RTCA ePacer station allows the cells to be continuously monitored. The E-Plate can also be inserted into the RTCA CardioECR station for monitoring.

Maintaining Cardiomyocytes in the E-Plate: Days 1 - 6

1. Equilibrate an aliquot of Maintenance Medium in a 37°C water bath immediately before use.
2. Maintain the cardiomyocytes on the E-Plate by replacing 100% of the spent medium with Maintenance Medium every 48 hours (100 µl/well).
 - a. Tilt the E-Plate to remove medium and gently add medium down the side of the well to avoid disturbing the cells.
3. Culture the cardiomyocytes in the RTCA ePacer at 37°C, 5% CO₂.

Pacing the Cardiomyocytes at 1.0 Hz: Day 7

Do not allow the pipettor tips to touch the bottom of the well during medium removal or addition. Tilt the E-Plate to remove medium and gently add medium down the side of the well to avoid disturbing the cells.

1. Equilibrate an aliquot of Maintenance Medium to 37°C.
2. Replace 100% of the spent medium with Maintenance Medium (100 µl/well).

- a. Tilt the E-Plate to remove medium and gently add medium down the side of the well to avoid disturbing the cells.
3. Insert E-Plate into the RTCA ePacer.
4. Initiate pacing by applying an electrical stimulus with a frequency of 1.0 Hz, pulse amplitude of 1.1 V and pulse width of 0.4 ms.

Note: *The electrical stimulus must contain a sufficient combination of sweeps and intervals to pace for at least 120 hours (e.g. 25 sweeps of 8 hour intervals).*

5. Maintain the cardiomyocytes on the E-Plate by replacing 100% of the spent medium with Maintenance Medium every 24 hours.
6. Culture the cardiomyocytes in the RTCA ePacer at 37°C, 5% CO₂ for 5 days.

Pacing the Cardiomyocytes at 1.5 Hz: Day 12

1. Equilibrate an aliquot of Maintenance Medium in a 37°C water bath.
2. Abort the previously applied electrical stimulus.
3. Replace 100% of the spent medium with Maintenance Medium (100 µl/well).
 - a. Tilt the E-Plate to remove medium and gently add medium down the side of the well to avoid disturbing the cells.
4. Insert E-Plate into the RTCA ePacer.
5. Initiate pacing by applying an electrical stimulus with a frequency of 1.5 Hz, pulse amplitude of 1.1 V and pulse width of 0.4 ms.

Note: *The electrical stimulus must contain a sufficient combination of sweeps and intervals to pace for at least 120 hours (e.g. 25 sweeps of 8 hour intervals).*

6. Maintain the cardiomyocytes on the E-Plate by replacing 100% of the spent medium with Maintenance Medium every 24 hours.
7. Culture the cardiomyocytes in the RTCA ePacer in a cell culture incubator at 37°C, 5% CO₂ for 5 days.

Pacing the Cardiomyocytes at 2 Hz: Day 17

1. Equilibrate an aliquot of Maintenance Medium in a 37°C water bath immediately before use.
2. Abort the previously applied electrical stimulus.
3. Replace 100% of the spent medium with Maintenance Medium (100 µl/well).
 - a. Tilt the E-Plate to remove medium and gently add medium down the side of the well to avoid disturbing the cells.
4. Insert E-Plate into the RTCA ePacer.
5. Initiate pacing by applying an electrical stimulus with a frequency of 2.0 Hz, pulse amplitude of 1.1 V and pulse width of 0.4 ms.

Note: *The electrical stimulus must contain a sufficient combination of sweeps and intervals to pace for at least 120 hours (e.g. 25 sweeps of 8 hour intervals).*

6. Maintain the cardiomyocytes on the E-Plate by replacing 100% of the spent medium with Maintenance Medium every 24 hours.
7. Culture the cardiomyocytes in the RTCA ePacer in a cell culture incubator at 37°C, 5% CO₂ for 5 days.

Measuring the Force-Frequency Relation of the Cardiomyocytes: Day 21+

1. Culture the cardiomyocytes in the RTCA ePacer at 37°C, 5% CO₂ without electrical stimulus for 4 hours to allow the culture to resume spontaneous beating rhythms.
2. Record a spontaneous beating measurement according to the RTCA ePacer Instrument Operator's Guide.
3. Pace the cells at 1.0 Hz for 1 sweep with a 10 minute interval. Enable the Auto function to perform the recording at the end of the step.
4. Pace the cells at 1.5 Hz for 1 sweep with a 10 minute interval. Enable the Auto function to perform the recording at the end of the step.
5. Pace the cells at 2.0 Hz for 1 sweep with a 10 minute interval. Enable the Auto function to perform the recording at the end of the step.
6. Evaluate the slope of the force-frequency relationship.
 - a. If the slope is positive, proceed to Data Acquisition and Analysis.
 - b. If the slope is negative, continue to apply electrical stimulus with a frequency of 2.0 Hz, pulse amplitude of 1.1 V and pulse width of 0.4 ms for approximately 24 hours before reevaluating the force-frequency relationship.




Figure 1: Representative Force-Frequency Relationship

iCell Cardiomyocytes², 01434 exhibit a negative force frequency relationship (black) at day 7 post-thaw. After pacing, the force-frequency relationship is positive (green) at day 22 post-thaw.

Data Acquisition and Analysis

Applying Compounds

1. Equilibrate an aliquot of Maintenance Medium to 37°C.

- Remove the spent medium and gently add 90 μl /well of Maintenance Medium to the side of the well to avoid disrupting the monolayer.
- Culture the cardiomyocytes in the RTCA CardioECR or ePacer at 37°C, 5% CO₂ without electrical stimulus to allow the culture to resume spontaneous beating rhythms.

Note: Spontaneous beating rhythms may not resume for up to 4 hours.

- Monitor the activity of the E-Plate to ensure a regular beating rhythm and stable whole-peak amplitudes are reached.
- Prepare test compounds in Maintenance Medium at 10X the final concentration in a 96-well cell culture plate.

Note: Final DMSO concentrations above 0.1% should be used with caution. Therefore the 10X compound solutions should not exceed 1% DMSO.

- Equilibrate the plate containing the 10X compound solutions at 37°C, 5% CO₂ for at least 1 hour.
- Quickly transfer 10 μl of the 10X compound solutions to the appropriate wells of the E-Plate.

Note: Beat rate and amplitude are temperature dependent. Do not keep the E-Plate outside the incubator for more than 5 minutes while compounds are added.

- Return the E-Plate to the RTCA CardioECR station. Record measurements according to the RTCA CardioECR Instrument Operator's Guide.

Data Acquisition and Analysis using the RTCA Software

See the RTCA CardioECR and ePacer Instrument Software Guides for specific instructions on using the software for data acquisition and analysis.

Example Data




Figure 2: Representative Effects of Long-term Electrical Pacing on Impedance

Paced iCell Cardiomyocytes², 01434 display the expected positive chronotropic response to Isoproterenol depicted by increased beat frequency and increased amplitude (green) when compared to a vehicle control (black).




Figure 3: Representative Effects of Electrical Pacing on Isoproterenol Response

iCell Cardiomyocytes², 01434 dose response curves for β -adrenergic agonist isoproterenol is concentration dependent. Panel A shows the expected positive chronotropic response when stimulated with isoproterenol. Panel B shows electrical pacing reverses the negative inotropic response to isoproterenol resulting in a greater than 3-fold increase in amplitude.

Customer's Responsibilities

FUJIFILM Cellular Dynamics, Inc. (FCDI), does not guarantee that you will obtain equivalent results from using iCell or MyCell products as described herein or that such use will not infringe any intellectual property right(s) of any third party(ies). You are solely responsible for obtaining any licenses you may require for your specific research use(s) of the iCell or MyCell products not expressly conveyed under FCDI's terms and conditions of sale or other transfer of the iCell or MyCell products to you.

Conditions of Use

The cells are FOR RESEARCH USE ONLY. See <https://fujifilmcdi.com/assets/tnc/standard.pdf> for USE RESTRICTIONS applicable to the cells and other terms and conditions related to the cells and their use.

Trademarks

iCell and MyCell are registered trademarks, and Cellular Dynamics and the  logo are trademarks of FUJIFILM Cellular Dynamics, Inc. All other brands, product names, company names, trademarks, and service marks are the properties of their respective owners.

Copyright Notice

© 2019 FUJIFILM Cellular Dynamics, Inc. All rights reserved. This document may not be reproduced, distributed, modified or publicly displayed without the express written permission of FUJIFILM Cellular Dynamics, Inc.

Revision History

Version 1.1: December 2019